COATEST® APC™ Resistance
Control Plasma Level 1 - 82 2650 63
 ENGLISH - Insert revision 01/2004

Intended use

For control of accuracy when using the Coatest APC Resistance and Coatest APC Resistance V kits.

Reagent

Control Plasma Level 1

Lympholized preparation of citrated human plasma.

CAUTION: Each donor unit used in the preparation of human source reagent has been tested by FDA approved methods for the presence of Hepatitis B surface antigen and antibodies to HIV 1 and 2 and Hepatitis C and found to be negative. However, since no test can completely rule out the presence of these blood borne diseases, the handling and disposal of human source reagents from this product should be made with care. Avoid contact with skin and eyes (S24/25). Do not empty into drains (S29). Wear suitable protective clothing (S36). This product is for *in vitro* diagnostic use.

Procedure

Reconstitute each vial with 1.0 mL deionized water, filtered through 0.22 µm or NCCLS type II water. Allow to stand at room temperature for 30 minutes. Swirl gently before use.

Storage conditions and stability

Stability after reconstitution is 6 hours at 2-25°C or 3 months at -20°C or below. The frozen reagent should be thawed at 37°C and gently mixed before use. Do not refreeze.

Expected values

Different clotting times will be obtained with different types of instruments depending on the clot detection principle. The following ranges cover the expected clotting times when performing analysis of Control Plasma Level 1 with the above batch number using the Coatest APC Resistance and Coatest APC Resistance V Kits.

Traceability of calibrators and control materials

The reported values were determined over multiple runs on ACL Futura using a specific lot of reagent and against an internal House Standard. As an international Standard is not available for the APCr-V assay, the values have been assigned against a House Standard which is traceable to frozen plasma pool.

Symbols used / Verwendete Symbole / Símbolos utilizados / Símbolos utilizados / Símboli impiegati / Símbolos utilizados / Anvendte symboler / Använda Symboler / Χρησιμοποιούμενα σύμβολα

ACCEPTANCE RANGE

Acceptance range / Referenzbereich / Rango esperado / Valeur usuelles / Intervallo di accettabilità / Acceptat område / Αποδεκτό εύρος

TP (e.g. AC, ML)

Turbidimetric/photometric / Turbidimetrisch/photometrisch / Turbidimétrico/fotométrico / Turbidimétrico/fotométrico / Turbidimétrico/fotométrico / PORTUGÜES / DANSK / SVENSK / GREEK

EM (e.g. SI, K, SI)

Electro-mechanical / Elektromechanisch / Electromecánico / Electromechanique / Elettro-mecanico / PORTUGÜES / DANSK / SVENSK / GREEK

ALL Instruments

All instruments / Alle Instrumente / Todos los instrumentos / Tous instruments / Tutti gli strumenti / PORTUGÜES / DANSK / SVENSK / GREEK

In vitro diagnostic medical device

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