

COAMATIC® Protein S-Free - 82 4003 63

Intended use

Automated latex ligand immunoassay for the quantitative determination of free Protein S (PS) in human citrated plasma on automated instruments.

Background and summary

Protein S is a vitamin K-dependent cofactor for the anticoagulant and the profibrinolytic effects of activated Protein C. Two forms of Protein S are present in plasma: free Protein S (40%), and Protein S linked to the complement C4b-binding protein (C4BP) (60%). Only free Protein S has functional cofactor activity.

Protein S deficiency may be hereditary³ or acquired. Acquired deficiency may be observed during pregnancy, oral anticoagulant therapy, oral contraceptive use, in liver disease, in newborn infants as well as in other clinical conditions. Deficiency of Protein S has been associated with a high risk of developing venous thromboembolism especially in young people.

Measurement principle

The COAMATIC® Protein S-Free assay determines the presence of free Protein S by measuring the increase of turbidity produced by the agglutination of two latex reagents. Purified C4BP adsorbed onto the first latex reagent, reacts with a high affinity for free Protein S of patient plasma in the presence of Ca²⁺ ions. The free Protein S adsorbed on the C4BP latex triggers the agglutination reaction with the second latex reagent, which is sensitized with a monoclonal antibody directed against human Protein S. The degree of agglutination will be directly proportional to the free Protein S concentration in the test sample.

Reagents

The COAMATIC® Protein S-Free kit consists of:

C4BP Buffer: 3 vials of Borax buffer containing bovine serum albumin, stabilizers and preservative.

C4BP Latex: 3 vials of a lyophilized suspension of polystyrene latex particles coated with purified human C4BP containing bovine serum albumin, stabilizers and preservative.

Anti PS MAb Latex: 3 vials of a suspension of polystyrene latex particles coated with a monoclonal antibody directed against human Protein S containing bovine serum albumin, stabilizers and preservative.

CAUTION:

The material in this product was tested by FDA approved test methods and found nonreactive for Hepatitis B Surface Antigen (HBsAg), Anti-HCV and HIV 1/2 antibodies. Handle as if potentially infectious.

PRECAUTIONS AND WARNINGS:

Avoid contact with skin and eyes (S24/25). Do not empty into drains (S29).

Wear suitable protective clothing (S36).

This product is for *in vitro* diagnostic use.

Reagent required but not provided

- | | |
|--------------------------------------|-------------------|
| 1. Calibration plasma | Art No 82 3534 63 |
| 2. Abnormal control plasma Level 1/2 | Art No 82 3559 63 |
| 3. Saline (0.9% NaCl) | |

Reagent preparation

C4BP Buffer: The reagent is ready for use.

C4BP Latex: Dissolve the contents of each vial of C4BP latex by pouring the entire contents of one vial of C4BP buffer into the latex reagent vial. Replace the stopper and swirl gently for a minimum of 20 seconds to completely dissolve the lyophilized latex. Make sure of the complete reconstitution of the product. It must appear as a homogenous and slightly milky suspension. Keep the reagent at 15-25°C for 30 minutes and invert to mix before use.

Do not shake.

Anti PS MAb Latex: Invert to mix before to use. Do not shake and AVOID FOAM FORMATION.

Storage condition and stability

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.

C4BP Latex - Stability after reconstitution: 1 month at 2-8°C in the original vial. Do not freeze.

Anti PS MAb Latex - Opened reagent is stable 1 month at 2-8°C. Do not freeze.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

Specimen collection and preparation

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Refer to NCCLS Document H21-A3 for further instructions on specimen collection, handling and storage.

Frozen plasma samples should be rapidly thawed at 37°C while gently mixing before testing. After thawing the assay must be performed within 2 hours.

Quality control

Normal and abnormal controls are recommended for a complete quality control program.¹³ Abnormal control plasma Level 1/2 Art No 82 3559 63 are designed for this program. The assigned values of these Controls are traceable to the International Standard. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance to good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard *et al* for identification and resolution of out-of-control situations.

Results

Free Protein S results are reported in % normality. The assay results should be used with other information, including the clinical context, in forming a diagnosis.

Limitations/interfering substances

Free Protein S results are not affected by heparin (UF heparin or LMW heparin) up to 1.5 IU/mL, bilirubin up to 18 mg/dL, hemoglobin up to 200 mg/dL, lipids up to 1280 mg/dL, platelets up to 10⁶ /L and rheumatoid factor up to 350 IU/mL. Hemolyzed and turbid samples should not be assayed.

COAMATIC® Protein S-Free assay is not affected by Factor V Leiden mutation (APC-R).

Reference values

A normal range study was performed using the COAMATIC® Protein S-Free kit on ACL Futura.

Sex	N	ACL Futura % fFree PS
Male	130	64.4 - 128.8
*Female	102	53.2 - 109.1

* Note: Age and hormonal status may affect the normal range for females.¹⁵

Ranges were calculated as recommended by the International Federation of Clinical Chemistry (IFCC). These results were obtained using a specific lot of reagent. Due to many variables, which may affect results, each laboratory should establish its own free Protein S normal range.

Calibration

A standard curve is obtained by analyzing different dilutions of Calibration Plasma Art. No 82 3534 63 in saline, which should be traceable to the International Standard.

ENGLISH - Insert revision 05/2003



Procedures

All conditions included in this package insert are referred to ACL Futura. Detailed instrument application setting including instructions for preparation of the reagents for a variety of automated instruments are available on request from Chromogenix.

Calculation

The changes in absorbance for the standards are plotted against free Protein S %. The samples results are then calculated from the linear equation obtained from the standard curve. This procedure is automatically handled by the instruments.

Performance characteristics

Precision:

Within run precision was assessed over multiple runs.

ACL Futura	CV% (Within run)	CV% (Between run)
Mean (%FPS)		
97.7 (n=60)	1.9	0.8
57.0 (n=60)	2.7	1.1
22.4 (n=60)	2.9	1.0

Correlation:

System	slope	intercept	r	Comparative method
ACL Futura	0.918	2.32	0.981	Free Protein S EIA
ACL 9000	0.933	5.09	0.972	Free Protein S EIA

The precision and correlation results were obtained using specific lots of reagents and controls.

Linearity: System

ACL Futura 12 - 135 (% Free Protein S)

Sample results above 135% should be manually diluted 1:2 with saline (200 µL of sample + 200 µL of saline) and reassayed. The printed results must be multiplied by 2 to correct for the dilution.

Detection Limity:

System

ACL Futura: 6.4 % Free Protein S

Determinations/Kit

On ACL Futura 75 tests (approx.)



Bibliography / Literatur / Bibliografía / Bibliographie / Bibliografia /Bibliografia / Litteratur / Litteraturförteckning / Βιβλιογραφία

1. Walker FJ. Regulation of Activated Protein C by a New Protein: A Possible Function for Bovine Protein S, J. Biol. Chem. 1980; 255: 5521-5524.

2. Suzuki K, Nishioka J. Plasma Protein S Activity measured using Protac, a Snake Venom Derived Activator of Protein C, Thromb. Res. 1988; 49: 241-251.

3. Sala N, Morell M, Tirado I, Espinosa Y, Llobet D, Fontcuberta J, Soria JM, Volpini V, Estivill X.. Linkage disequilibrium between the protein S Heerlen allele and protein S deficiency in Spanish families. Thromb. Haemost. 1995; 73: 1259.

4. Mannucci PM, Valsecchi C, Krachmalnicoff A, Faioni EM, Tripodi A. Familial Dysfunction of Protein S, Thromb. Haemost. 1989; 62: 763-766.

5. D'Angelo A, Viganò D'Angelo S, Esmon CT, Comp PC. Acquired Deficiencies of ProteinS. Protein S Activity during Oral Anticoagulation, in Liver Disease, and in Disseminated Intravascular Coagulation, J. Clin. Invest. 1988; 81: 1445-1454.

6. Comp PC, Esmon CT. Recurrent Venous Thromboembolism in Patients with a Partial Deficiency of Protein S, N. Engl. J. Med. 1984; 311: 15252-1528.

7. Maccaferri M, Legnani C, Preda L, Palareti G. Protein S Activity in Patients with Heredofamilial Protein S Deficiency and in Patients with Juvenile Venous Thrombosis. Results of a functional method, Thromb. Res. 1991; 64: 647-658.

8. Faioni EM, Valsecchi C, Palla A, Taioli E, Razzari C, Mannucci PM. Free Protein S Deficiency is a Risk Factor for Venous Thrombosis. Thromb. Haemost. 1997; 78: 1343-1346.

9. Dahlbäck B. Purification of human C4b-binding protein and formation of its complex with vitamin K-dependent protein S, Biochem. J. 1983; 209: 847-856.

10. He XH, Shen L, Malmborg AC, Smith KJ, Dahlbäck B, Linse S. Binding site for C4b-binding protein in vitamin K-Dependent protein S fully contained in carboxy-terminal laminin-G-type repeats. A study using recombinant factor IX protein S chimeras and surface plasmon resonance. Biochemistry 1997; 36: 3745-3754.

11. Richmond JY, McKinney RW eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, 4th Edition, 1999.






12. National Committee for Clinical Laboratory Standards. Collection, Transport and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays, Third Edition, NCCLS Document H21-A3; Vol. 18 No. 20.

13. Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, Am. J. Clin. Pathol. 1970; 53: 924-927.

14. Westgard JO, Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Processes, AACC Press 1986.

15. Liberti G, Bertina RM, Rosendaal FR. Hormonal State rather than Age Influences Cut-off Values of Protein S: Reevaluation of the Thrombotic Risk Associated with Protein S Deficiency, Thromb. Haemost. 1999; 82: 1093-1096. Solberg H.E. Approved Recommendation (1987) on the Theory of Reference Values. Part 5. Statistical Treatment of Collected Reference Values. Determination of Reference Limits. J. Clin. Chem. Clin. Biochem. 1987; 25: 645-656.

Symbols used / Verwendete Symbole / Símbolos utilizados / Symboles utilisés / Simboli impiegati / Símbolos utilizados / Anvendte symboler / Använda Symboler / Χρησιμοποιηθέντα σύμβολα

<div>IVD</div> <div><div>In vitro diagnostic medical device</div><div>In-vitro Diagnostikum</div><div>De uso diagnóstico <i>in vitro</i></div><div>Dispositif médical de diagnostic <i>in vitro</i></div><div>Per uso diagnostico <i>in vitro</i></div><div>Dispositivo médico para utilização em diagnóstico <i>in vitro</i></div><div>"in vitro" diagnostisk udstyr</div><div><i>In vitro</i> diagnostisk medicinsk produkt</div><div>Προϊόν για διαγνωστική χρήση <i>In vitro</i></div></div>	<div>LOT</div> <div><div>Batch code</div><div>Chargen-Bezeichnung</div><div>Identificación número de lote</div><div>Désignation du lot</div><div>Numero del lotto</div><div>Número de lote</div><div>Batch nr.</div><div>Tillverkningskod</div><div>Αρ. Παρτίδας</div></div>	<div></div> <div><div>Use by</div><div>Verwendbar bis</div><div>Caducidad</div><div>Utilisable jusqu'à</div><div>Da utilizzare prima del</div><div>Data limite de utilização</div><div>Anvendelse</div><div>Användning</div><div>Χρήση έως</div></div>	<div></div> <div><div>Temperature limitation</div><div>Festgelegte Temperatur</div><div>Temperatura de Almacenamiento</div><div>Températures limites de conservation</div><div>Limiti di temperatura</div><div>Limite de temperatura</div><div>Temperatur begrænsninger</div><div>Temperatur gräns</div><div>Περιορισμοί θερμοκρασίας</div></div>	<div></div> <div><div>Consult instructions for use</div><div>Beilage beachten</div><div>Consultar la metódica</div><div>Lire le mode d'emploi</div><div>Vedere istruzioni per l'uso</div><div>Consultar as instruções de utilização</div><div>Se vejledning for anvendelse</div><div>Ta del av instruktionerna före användning</div><div>Συμβουλευτήτε τις οδηγίες χρήσης</div></div>	<div>CONTROL</div> <div><div>Control</div><div>Kontrollen</div><div>Control</div><div>Contrôle</div><div>Controllo</div><div>Controlo</div><div>Kontrol</div><div>Kontroll</div><div>Υλικό ποιοτικού ελέγχου</div></div>	<div></div> <div><div>Biological risks</div><div>Biologisches Risiko</div><div>Riesgo biológico</div><div>Risque biologique</div><div>Rischio biologico</div><div>Risco biológico</div><div>Miljø oplysninger</div><div>Biologiska risker</div><div>Βιολογικοί κίνδυνοι</div></div>	<div></div> <div><div>Manufacturer</div><div>Hergestellt von</div><div>Fabricado por</div><div>Fabricant</div><div>Prodotto da</div><div>Fabricado por</div><div>Producent</div><div>Tillverkare</div><div>Κατασκευαστής</div></div>	<div><div>EC</div><div>REP</div></div> <div><div>Authorised representative</div><div>Bevollmächtigter</div><div>Representante autorizado</div><div>Mandataire</div><div>Rappresentanza autorizzata</div><div>Representante autorizado</div><div>Leverandør</div><div>Auktoriserad representant</div><div>Εξουσιοδοτημένος αντιπρόσωπος</div></div>
--	--	---	--	--	--	--	---	--