



ELISA Fibronectin

Complete ELISA kit for the quantitative determination of human fibronectin.

Test principle

The samples are incubated in micro wells coated with an anti-fibronectin antibody. After binding of fibronectin to the capture antibody, unbound substances are removed by washing. A peroxidase conjugated detecting antibody is added and binds to fibronectin /capture antibody complex. After a washing step, added substrate reacts with the conjugate. The intensity of the resulting colour is proportional to the content of fibronectin in the sample.

Reagents and stability (after opening)

Micro wells	12x8	8 weeks at 4°C
Diluent buffer	3 vials	4-8 weeks at 4°C
Calibrator	1 vial	4 hours at RT 8 weeks at -70°C
Conjugate	1 vial	4 hours at RT
Wash buffer (12.5 x)	1 vial	4-8 weeks at 4°C
TMB-substrate	1 vial	until expiry date at 2-8°C
Stop solution	1 vial	until expiry date at 2-8°C

No. of determinations: 96

Article no.: TC12030