



Emicizumab Controls Level 1 and Level 2

Emicizumab Level 1 and Level 2 controls are used for quality control of modified one-step FVII activity assays to determine the concentration of active emicizumab from a clotting time measurement.

Test principle

Level 1 and Level 2 emicizumab controls can be used with all tests as with citrated plasma samples. They are used to monitor test variables for the one-step FVIII activity assay to determine the concentration of active emicizumab.

Reagent and stability (reconstituted)

The reagent must be stored unopened at +2-8°C and can be used until the date indicated on the label.

Emicizumab Control Level 1	5 x 1 ml	8 h at 2-8 °C 8 h at 23-25 °C
Emicizumab Control Level 2	5 x 1 ml	8 h at 2-8 °C 8 h at 23-25 °C

Article No.:

152-401-CE