






UNITEST Factor XI

A chromogenic method for the determination of factor XI activity in plasma.

Test principle



Reagents and stability (after opening)

Factor XIa substrate	1 vial	6 months at 4°C
Kaolin	1 vial	see expiry date
Buffer concentrate	1 vial	see expiry date
Kallikrein inhibitor	1 vial	8 hours at 4°C
Corn trypsin inhibitor	1 vial	8 hours at 4°C
Standard plasma	1 vial	8 hours at 4°C

Article no.: 0090