



Technoplastin® HIS

Coagulation test for determination of thromboplastin time. The lyophilized reagent from rabbit brain contains a heparin neutralizer. ISI ~ 1,2

Test principle

The sample warmed to 37°C is mixed with Technoplastin HIS. The time until formation of a clot is measured.

Reagents and stability (after opening)

Technoplastin HIS

- 6 hours at 37°C
- 24 hours at RT
- 6 days at 2-8°C
- 2 months at -20°C

Analyzer applications:

on request

Package size: 12 x 2 ml

Article no.: 5003009

Package size: 2 x 10 ml

Article no.: 5003030

Package size: 6 x 10 ml

Article no.: 5003026

Package size: 20 x 10 ml

Article no.: 5003021