



UNITEST Factor XII

A chromogenic method for determination of the factor XII-level in plasma.

Test principle

Factor XII $\xrightarrow{\text{Activator}}$ Factor XIIa

Substrate $\xrightarrow{\text{Factor XIIa}}$ Peptide + pNA

Reagents and stability (after opening)

Factor XII-Substrate	1 vial	6 months at 4°C
Factor XII-Activator	1 vial	8 hours at 4°C
		6 months at <-20°C
Buffer concentrate	1 vial	s. expiry date
Kallikrein-Inhibitor	1 vial	8 hours at 4°C
		6 months at <-20°C
Standard plasma	1 vial	8 hours at 4°C

Article no.: 0068