



Technoclot® PT Plus

Coagulation test for determination of thromboplastin time. The lyophilized reagent from rabbit brain contains a heparin neutralizer. ISI ~ 1,1

Test principle

The sample warmed to 37°C is mixed with Technoclot PT. The time until formation of a clot is measured.

Reagents and stability (after opening)

Technoclot PT	6 hours at 37°C
	24 hours at RT
	6 days at 2-8°C
	2 months at -20°C

Analyzer applications:

on request

Package size:	12 x 4 ml
---------------	-----------

Article no.:	5003204
--------------	---------

Package size:	20 x 10 ml
---------------	------------

Article no.:	5003210
--------------	---------