



UNITEST Kallikrein-like Activity

A chromogenic method for the determination of kallikrein-like activity in plasma, which is mainly due to the activity of α_2 Macroglobulin-bound kallikrein.

Test principle

Substrate $\xrightarrow{\alpha_2\text{M} - \text{Kallikrein}}$ Peptide + pNA

Reagents and stability (after opening)

Kallikrein-Substrate	1 vial	6 months at 2-8°C
Buffer-Concentrate	1 vial	see expiry date
Standard plasma (high activity)	1 vial	4 hours at RT
Standard plasma (low activity)	1 vial	4 hours at RT

Article no.:

0087