




## UNITEST Kallikrein-Inhibitors

A chromogenic method for the determination of kallikrein inhibitors in plasma.

### Test principle

Kallikrein-INH + Kallikrein (excess) →

[Kallikrein-INH • Kallikrein] + Kallikrein (residual)



### Reagents and stability (after opening)

Kallikrein-Substrate	1 vial	8 hours at 2-8°C 6 months at -20°C
Kallikrein (human)	1 vial	see expiry date
Buffer-concentrate	1 vial	see expiry date
Standard plasma	1 vial	8 hours at 2-8°C

Article no.:

**0073**