



Thrombin Reagent

Coagulation test for determination of thrombin time. The lyophilized reagent contains per vial 20 international units bovine thrombin.

Test principle

The sample warmed to 37°C is mixed with the thrombin reagent. The time until formation of a clot is measured.

Reagents and stability (after opening)

Thrombin Reagent	8 hours at RT
	1 Month at -20°C

Package size:	6 x 6 ml
----------------------	-----------------

Article no.:	5100005
---------------------	----------------

Analyzer applications:	on request
-------------------------------	-------------------