



Endotoxin-free related products

Haemotox rFC related endotoxin-free products.

Haemotox Tubes

Borosilicate-glass test tubes with aluminium screw cap. Applied to dilution of endotoxin standard as well as to dilution and aliquotation of samples. LPS content: < 0,005 EU/ml

Capacity	5 ml
Height	50 mm
Outer diameter	16,1 mm
Thickness	0,95 mm
Amount	112

Article no.: **31221**

Haemotox Water

Sterile, ultra-pure water. Recommended for endotoxin-free reconstitution and dilutions of reagents, samples and standards used in Endotoxin analytics. LPS content: < 0,001 EU/ml

Haemotox Water 30 ml 20 x 30 ml

Article no.: **31222**

Haemotox Water 100 ml 16 x 100 ml

Article no.: **31223**

Haemotox Water 500 ml 6 x 500 ml

Article no.: **31224**

Endotoxin-free containers

Certified endotoxin-free containers for sample storage.

30 ml, unlabeled, ca. 400 containers

Article no.: **AB128**

30 ml, labeled, ca. 400 containers

Article no.: **AB128L**

100 ml, unlabeled, ca. 200 containers

Article no.: **AP185**

100 ml, labeled, ca. 200 containers

Article no.: **AP185L**