



COAMATIC® Protein C

A chromogenic method for the determination of Protein C in plasma. Not influenced by heparin levels up to 3 IU/ml.

Test principle

Protein C $\xrightarrow{\text{Protein C Activator}}$ APC

S-2366 $\xrightarrow{\text{APC}}$ Peptide + pNA

Reagents and stability (after opening)

S-2366	2 vials	3 months at 2-8°C
Protein C Activator	2 vials	3 months at 2-8°C

Analyzer applications:

On request

Article no.:

82 20 98