



Heparin Resistant Calcium Chloride

For discrimination of heparin from oral anticoagulants, factor deficiencies and lupus coagulants in both PT and aPTT systems.

Test principle

Addition of a heparin neutralizing agent to the calcium chloride solution rather than the sample or reagent avoids interference with contact activation and allows discrimination of heparin from other anticoagulants, factor deficiencies and lupus anticoagulants.

Reagent and stability (after opening)

Recalcifying Solution

expiry date

Package size:

10 x 5 ml

Article no.:

X9107-HK