



COAMATIC® Factor VIII

Chromogenic method for determination of factor VIII activity in plasma, blood fractions and concentrates. Fulfils the requirements of the European Pharmacopoeia for factor VIII concentrate testing.

Test principle

Factor X $\xrightarrow{\text{FIXa, CaCl}_2, \text{phospholipids, FVIII}}$ Factor Xa

S-2765 $\xrightarrow{\text{Factor Xa}}$ Peptide + pNA

Reagents and stability (after opening)

S-2765 + I-2581	1 vial	1 month at 2-8°C
Factor reagent	2 vials	1 day at 2-8°C 1 month -70°C
Buffer	1 vial	1 month at 2-8°C

Analyzer applications:

on request

Article no.:

82 25 85