



LPC Software

Software for controlling APACT 4004

Simple and flexible for routine use

- User friendly and flexible software with graphical menus
- Free and adaptable presentation of reaction curves from all connected measuring channels
- Possibility to connect 1 to 6 APACT 4004, enabling up to 24 simultaneous channels with USB connections
- Data are automatically saved in a date organized database
- Coloured presentations of inductors and curves, with independent modification options
- Task presentation, assay names and concentrations are editable
- Performance of ristocetin-cofactor-assays (optional)
- Automatic evaluation of max. and min. aggregation, shape change, lag phase, area under the curve, disaggregation and further parameters
- Result calculation and evaluation options free for editing
- Multi-language software

Article no.:

30.000.3778